

LONG-TERM CARE BEDS




WE CARE LIKE FAMILY™

Reference Chart




Unique Features	Bed Type					
	12001: D100	12002: D200	12003: D300	12004: D400	12005: D500	12013: DB300
Split Frame	—	—	—	Yes	Yes	—
Compares to Drive Medical	None	P503	P703 - P903	N/A	N/A	P750
Bed Height	10.5" - 25"	7.75" - 30"	7" - 30"	4" - 22"	2" - 22"	10.5" - 30"
Bed Width	36"	36"	36" - 42"	36"	36"	36" - 42" - 48"
Width Extender	—	—	Optional	—	—	Included
Deck Length	80"	80"-84"-88"	80"-84"-88"	80"	80"	80"-84"-88"
Embedded Length Extender	—	Yes	Yes	—	—	Yes
Quick Release Boards 36"	Yes	Yes	Yes	—	—	Yes
Quick Release Adjustable Width Headboard (Optional)	—	—	36" - 42"	—	—	36" - 48"
DC Motors	3	3	4	4	4	4
Hand Control Functions	4 - Head, Foot, Hi-Lo, Auto Contour	4 - Head, Foot, Hi-Lo, Auto Contour	5 - Head, Foot, Hi-Lo, Auto Contour, Trend and Reverse Trend	5 - Head, Foot, Hi-Lo, Auto Contour, Trend and Reverse Trend	5 - Head, Foot, Hi-Lo, Auto Contour, Trend and Reverse Trend	5 - Head, Foot, Hi-Lo, Auto Contour, Trend and Reverse Trend
Auto Contour	Yes	Yes	Yes	Yes	Yes	Yes
Trendelenburg/Reverse Trendelenburg	—	—	Yes	Yes	Yes	Yes
Comfort Chair Position	—	—	Yes	Yes	Yes	Yes
Rail Type: Half, Assist Bar, Edge Bar	All	All	All	All	All	All
Wall Bumper Guard	Yes	Yes	Yes	—	—	Yes
Wheel Type and Central Locking System	4 Locking	4 Locking 4 Non-Locking	2 Locking at head and CLS at foot	4 Locking	4 Locking	2 Locking at head and CLS at foot
Optional Battery Back Up	Yes	Yes	Yes	Yes	Yes	Yes
Weight capacity	450 lb.	450 lb.	500 lb.	450 lb.	450 lb.	750 lb.




MH – Mahogany



CH – Amber Cherry



OA – Light Oak



CM – Composite
(for D200, D300 & DB300)