NitriDerm<sup>®</sup> CO

## Nitrile Exam Gloves Series 125

COATS® barrier technology is a totally powder free, all-natural physical barrier between the glove surface and the wearer's skin, utilizing the anti-irritant and anti-inflammatory properties of colloidal oatmeal. In our COATS® gloves, the colloidal oatmeal barrier not only shields the skin surface, it also creates a moisturizing and pH balancing environment that helps dry, damaged skin to repair itself. And healthy skin is the body's best defense against infection. What a naturally great idea!

Characteristics: Non-Sterile, Powder-Free

	Will the supplementary
- ORD	The state of the s
1 Comments	many improved view immers which
1 Comments	many improved view immers which
1 Comments	Tri Dermi COATS  Tree Nimit Lam Glavo Callidad Gunned Adire Therapeatic System  Grand for use will be compared to the control of the control
N. Powde	TriDerni COATS  TriDerni CoATS  TriDerni CoATS  TriDerni CoATS  TriDerni CoATS  Tride Territorial System  Tributorian General Control Heavy College Control Co
N. Powde	Fri Dermi COATS Free Nutrile Euron Glaves - Calabella Osmoral Joseph Chameled System  Special Control of Contr
N. Powde	Tree Name Exam Gloves . Collisida Outmood Lettive Thorappeatic System Free Name Exam Gloves . Collisida Outmood Lettive Thorappeatic System Free Name Exam Gloves . 125 mars G
1 Comments	Tri Derni COATS  Free Nimit Eum Glave Collidad Gunned Adive Therapeatic System  Committee Commit

Tested for use with Chemotherapy Drugs

SIZE	ITEM NO.	
XS	125052	
S	125102	
М	125202	
L	125302	
XL	125352	

## All Sizes: 200 Gloves/Box, 10 Boxes/Case

## Features:

- Flat Stacking = Gloves Dispense Easily, One by One
- Colloidal Oatmeal = Therapeutic
- Non-Latex = No Latex Allergens
- Finger Textured = An Improved Wet/Dry Grip
- Low Modulus = A Softer, More Comfortable Fit

Manufactured in compliance with multiple international standards, including ASTM D6319, ASTM D5151, ASTM F1671, and ASTM D6978

Tested for use with Chemotherapy Drugs as per ASTM D6978 Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs Minimum Breakthrough Detection Time in Minutes:						
Carmustine (3.3 mg/ml)	17.2	Fluorouracil (50.0 mg/ml)	>240			
Cisplatin (1.0 mg/ml)	>240	Methotrexate (25.0 mg/ml)	>240			
Cyclophosphamide (20.0 mg/ml)	>240	Mitomycin C (0.5 mg/ml)	>240			
Dacarbazine (10.0 mg/ml)	>240	Paclitaxel (6.0 mg/ml)	>240			
Doxorubicin Hydrochloride (2.0 mg/ml)	>240	Thiotepa (10.0 mg/ml)	36.1			
Etoposide (20.0 mg/ml)	>240	Vincristine Sulfate (1.0 mg/ml)	>240			

Do not use with Carmustine. Do not use with Thiotepa.

Caution: In compliance with ASTM D6978, the testing conditions used are intended to approximate the worst-case conditions for clinical use. Testing was conducted on single layer glove material. It is the users' responsibility to determine the applicability of these gloves for their intended use with chemotherapy drugs.

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength (Mpa)	34	32.6
ASTM Requirement Min. (Mpa)	14	14
Elongation (%)	632	583
ASTM Requirement Min. (%)	500	400

Average Length	9.5 in <b>◆</b> 240 mm
Average Palm Thickness	3.2 mil ◆ 0.08 mm
Average Finger Thickness	3.8 mil ◆ 0.10 mm





Product Solutions You Trust

View this product on our website:

