

Flat Top, Classic Series, Straight Line Treatment Table with Drawers





- **MODEL 1300** • Solid hardwood legs
- H-brace construction
- Two drawer unit and shelf in Natural or Dark Cherry laminate to coordinate with frame
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in Natural and Dark Cherry finish
- Choice of upholstery colors
- Fast, easy assembly with basic tools

Strap

• 400 lbs. load capacity under normal use (181.43 kg)

& Cutter Combo

SPECIFICATIONS

1300-27	72" L x 27" W x 31" H
1300-30	72" L x 30" W x 31" H
1300-38	78" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

Options



Finishes



Upholstery Colors

