

# DermAssist® Prestige®

## Latex Surgical Gloves Microsurgical Series 137

These light brown powder-free latex surgical gloves are 30% thinner than standard surgical gloves. They are ideal for use when maximum tactile sensitivity is required. The strong, yet highly sensitive latex barrier is roughened for improved control and sensitivity. They can also be used as the inner glove in a double-gloving system.

Characteristics: Powder-Free, Sterile



SIZE	ITEM NO.
5.5	137550
6	137600
6.5	137650
7	137700
7.5	137750
8	137800
8.5	137850
9	137900



### Features:

- Polymer Coated = Damp-hand Donnability
- Bisque Finish = an Improved Wet/Dry Grip
- Low Protein = Reduced Chances of Latex Allergic Reactions

All Sizes:  
25 Pairs/Box, 4 Boxes/Case

Manufactured in compliance with multiple international standards, including ASTM D3577, ASTM D5151, ASTM F1671, and D5712

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength (Mpa)	29	27
ASTM Requirement Min. (Mpa)	24	18
Elongation (%)	860	820
ASTM Requirement Min. (%)	750	560

Average Length	11.6 in ♦ 295 mm
Average Palm Thickness	5.9 mil ♦ 0.15 mm
Average Finger Thickness	6.3 mil ♦ 0.16 mm



**Intertek**

Innovative Healthcare Corporation is certified to ISO 13485:2016 QMS for medical devices.



**INNOVATIVE  
HEALTHCARE  
CORPORATION**

www.ihc-solutions.com  
1-800-272-1533

*Product Solutions You Trust*

View this product  
on our website:

