

DermAssist® Prestige®


Latex Surgical Gloves

DHD

Series 139

Damp don gel interior provides effortless donning. Curved finger design increases comfort. Powder-Free to reduce risk of powder induced allergies, adhesion, and granulomas.

Characteristics: Powder-Free, Sterile



SIZE	ITEM NO.
5.5	139550
6	139600
6.5	139650
7	139700
7.5	139750
8	139800
8.5	139850
9	139900

All Sizes:
50 Pairs/Box, 4 Boxes/Case



Features:

- Unique Hydrogel Interior = Damp-hand Donnability
- Bisque Finish = an Improved Wet/Dry Grip
- Low Protein = Reduced Chances of Latex Allergic Reactions

Manufactured in compliance with multiple international standards, including ASTM D3577, ASTM D5151, ASTM F1671, and D5712

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength (Mpa)	29	27
ASTM Requirement Min. (Mpa)	24	18
Elongation (%)	860	820
ASTM Requirement Min. (%)	750	560

Average Length	11.6 in ♦ 295 mm
Average Palm Thickness	7.5 mil ♦ 0.19 mm
Average Finger Thickness	8.7 mil ♦ 0.22 mm



Intertek

Innovative Healthcare Corporation is certified to ISO 13485:2016 QMS for medical devices.



**INNOVATIVE
HEALTHCARE
CORPORATION**

www.ihc-solutions.com
1-800-272-1533

View this product
on our website:



Product Solutions You Trust