

# Part Number/Price

Part #	Model #	Capacity	Price
121304	150-10-7	550 × 0.2 lb (250 × 0.1 kg) with USB	Consult
194728	150-10-7BLE	550 × 0.2 lb (250 × 0.1 kg) with Bluetooth® BLE 4.0	Consult

# **Shipping Weight**

Part #	Model #	L×W×H	Package Wgt. (lb)	Dim Wgt. (lb)	Shipping Wgt. (lb)
121304	150-10-7	25 × 18 × 6 in	17	20	20
194728	150-10-7BLE	25 × 18 × 6 in	17	20	20

# **Options/Accessories**

Part #	Description	Price
102508 ††	9 VDC adaptor, 230 V	Consult
107445	Carrying case	Consult
115655 <sup>†</sup>	9 VDC adaptor, 120 V/230 V auto switching c-MET medical grade power supply	Consult
133077 <sup>†</sup>	9 VDC adaptor, 120 V/230V auto switching	Consult
81076	Epson TM-U220 tape printer	Consult
186830	Printer cable, 25-pin D-Sub male to RJ-45, 6 ft	Consult






† NEMA 1-15 Plug

†† CEE 7/16 Plug

# **Dimensions**

	Inches	Centimeters
Α	14.7	37.5
В	14.7	37.5
С	3.0	7.6
D	7.8	20.0
E	6.7	17.0



# **Standard Features**

- · Large, stable scale base
- · Push-button LCD display
- Unit of measure toggle or locking
- Motiontrap™ movement compensation technology
- · EHR capable
- USB or Bluetooth® BLE 4.0 available, part number selectable
- · Hold function
- · Zero function
- BMI (body mass index) function
- 6 ft cable between digital display and base
- · Factory calibrated

# **Specifications**

**Capacity:** 550 × 0.2 lb (250 × 0.1 kg)

# **Platform Dimensions:**

(W × L × H) 14.5 × 14.5 × 3 in

## Power:

9 VDC, provided by six AA alkaline batteries (included) or AC adaptor (optional)

# **Operating Temperature:** 50 °F to 104 °F (10 °C to 40 °C)

# **Battery Type:**

Six AA alkaline (included)

## **Battery Use:**

25 hours continuous use with batteries Automatic power-off can be configured

5-digit LCD display, 0.75 in (1.9 cm) digit height

# Approvals:

c-Met, E113986

# Warranty:

Two-year limited