

Privacy Screens



#153094 w/Apricot - #153020

Economy Folding Screen Frame

#153094

The 4 sections of the screen are attached to each other by a riveted hinge which allows the frame to smoothly fold and unfold as needed. Provides instant isolation for routine, temporary and special circumstances. Anodized aluminum frame and "accordion" hinge allows frame to fold for convenient storage. Screen panels are secured to the frame by top and bottom snap in rods.

Features:

- Available in 2, 3 or 4 section models
- Four sections measures 68"H x 54"W
- Requires one vinyl or fabric panel per section
- Anodized aluminum construction provides a lightweight and durable solution eliminating the possibility of paint chipping or rusting
- 5 year warranty



White
#153010



Green
#153015

Screen Frames and Panels are Sold separately



Apricot
#153020




Palmetto
#153027



Oatmeal
#153028



Marine
#153029




Rainbow
#153031




Norway
#153035



Rouge
#153040




Autumn
#153042




Powder
#153044



Frost
#153045



Lavender
#153046



Sage
#153048