



Products / Personal Protection / Face Masks / Masks / Anti-Fog Surgical Mask with Tie 3-Ply, Blue



Anti-Fog Surgical Mask with Tie 3-Ply, Blue

1550

Tie, 50/BX, 6 BX/CS

Product Features

3-Ply pleated lightweight, breathable material

BFE of $\geq 99\%$

Meets ASTM Level 1 standard

Attributes

| | |
|----------------------|---------|
| Brand Name | Dukal |
| Color | Blue |
| HCPC Code | A4928 |
| Sterile | No |
| Weight | 4.4 lbs |
| Everyday Care | No |
| Premium Care | No |

Packaging

| | |
|--------------------------------|--------------------|
| Packaging | 50/BX, 6 BX/CS |
| Inner Packaging | 50/BX |
| Case Packaging | 6 BX/CS |
| Total Case Quantity | 300 |
| Item Packaging Barcode | |
| Inner Packaging Barcode | (01)10665973004980 |
| Case Packaging Barcode | (01)20665973004987 |

2 Fleetwood Ct.
Ronkonkoma, NY 11779

Phone: 1 (631) 656-3800
Toll Free: 1 (800) 243-0741
Fax: 1 (631) 656-3810

