

## Nitrile Powder Free Examination Glove - Chemo Rated

Our OmniTrust 212 Series Nitrile Powder Free Examination Gloves - Chemo Rated is a thermally activated glove that starts conforming to your hand when donned. You will find our Comfort Formulation very comfortable, reducing hand fatigue. It has our famous (micro-texture) fingertips for added grip control without loss of tactile sensitivity. Our #212 series is also tested safe use with Chemotherapy Drugs.



- Ambidextrous Beaded Cuff, Single Use Only
- Comfortable Formulation, Thermally Activated
- Tested For Safe Use With Chemo Drugs
- Not Made With Natural Rubber Latex
- Textured Fingers for Tactile Sensitivity
- Fentanyl Resistant

### Tensile Strength - Aged/Unaged

min. 14 MPa / min. 15 MPa

### Elongation - Aged/Unaged

min. 400% / min. 500%

Dimensions	Finger	Palm	Cuff	Length	Width
<b>X-Small</b>	min. 0.08 mm 3.15 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	75 ± 5 mm
<b>Small</b>	min. 0.08 mm 3.15 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	85 ± 5 mm
<b>Medium</b>	min. 0.08 mm 3.15 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	95 ± 5 mm
<b>Large</b>	min. 0.08 mm 3.15 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	106 ± 5 mm
<b>Extra Large</b>	min. 0.08 mm 3.15 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	116 ± 5 mm

Item Numbers	X-Small	Small	Medium	Large	X-Large
<b>212 Series</b>	212-20	212-21	212-22	212-23	212-24

Packaging	Boxes	Cases
packaged by weight	100 Gloves	10 Boxes

