


## TRAVELITE TRANSPORT CHAIR (TC005GY)



- |   |   |
|---|---|
| 1. Handle Bar Left- TC00501L<br>Handle Bar Right- TC00501                                 | 11. Footrest Assembly (pr)- TC005012                        |
| 2. Handgrip- TC00502  | 12. Rear Wheel Axle & Spacer- TC005014                      |
| 3. Upholstery- TC00503  | 13. Caster Axle- TC005015                                   |
| 4. Arm Assembly Left- TC00504L<br>Arm Assembly Right- TC00504R<br>Armrest Tube- TC005013L | 14. Cross Frame- TC005017                                   |
| 4a. Armrest Latch- TC00507  | 15. Front Frame Tube Right- TC005018                        |
| 4b. Arm Pad 10" Pad- TC00508<br>Arm Pad 8" Pad- TC005028                                  | 16. Front Frame Tube Left- TC005019                         |
| 5. Folding Hinge Left- TC00505L<br>Folding Hinge Right- TC00505                           | 17. Fork- TC005020<br>Fork Bearing- STDS1060                |
| 6. Telescoping Tube- TC00506  | 18. Rear Wheel- TC005021R                                   |
| 7. Cross Frame Block (4 ea)- TC005011   | 18a. Rear Wheel Bearing (pr)- STDS4S3602                    |
| 8. Upholstery Strap- TC005027   | 19. Caster- TC005021F                                       |
| 9. Footrest Connector- TC00509  | 19a. Caster Bearing- STDS1011                               |
| 10. Footrest Tube- TC005010   | 20. Wheel Lock Left- TC005022<br>Wheel Lock Right- TC005023 |
|   | 21. Tipping Sleeve- STDS3J2422                              |
|   | 22. Frame Plug- TC005026                                    |
|   | 23. Carry Bag- TC005024                                     |