Telescopic measuring rod with large measuring range



seca 222

- + Designed specifically for wall mounting
- + With side print for greater ease of reading
- + Extra wide head gauge
- + Includes heel positioner for ideal three-point measurement
- + Durable aluminum measuring rod
- + Abrasion-resistant scale due to high-quality screen printing



Telescopic measuring rod with large measuring range





Mechanical telescopic measuring rod with large measuring range

The seca 222 is specially designed for wall mounting. This classic seca personal measuring rod with a large measuring range of 2 to 90 inch and an additional heel positioner is selected by thousands of health care professionals worldwide due to its accuracy and intuitive operation.

An excellent measure

Thanks to the seca quality tested in continuous load tests, the seca 222 pull-out measuring rod will still deliver the most precise results even after years of daily use. The high-quality finish of the aluminum profile ensures that handling is particularly smooth. The printed scale is abrasion-resistant and the large head positioner is interchangeable and can be folded away when not in use.

Measuring range	2 – 90 in / 6 – 230 cm
Graduation	1/8 in / 1 mm
Dimensions (WxHxD)	12.1 x 91.3 x 11.3 in / 308 x 2,320 x 286 mm
Dimensions when folded (WxHxD)	12.1 x 52.2 x 1.5 in / 308 x 1,325 x 37 mm
Net weight	2.4 lbs / 1.1 kg



Height measurement with millimeter precision

Wall mounting makes the seca 222 a particularly stable and spacesaving measuring instrument. The foot positioner enables precise three-point measurement. Together with the reading marking at abdominal height and the abrasion-resistant scale, height is always measured with millimeter precision.