Non-Latex • Powder-Free • Textured















SemperCare® Nitrile combines superb barrier protection, comfort, fit, and feel while providing the tensile strength necessary for demanding environments. An easy to don glove, SemperCare® is an excellent latex-free option with textured fingertips and palm ensuring a safe grip, even in a wet environment. With an instantly recognizable blue nitrile color, this glove is an ideal choice for your entire team. At Sempermed®, we care about your comfort and protection. Always.

## **FEATURES**

- Excellent strength and barrier/chemical resistance
- Better elasticity
- Not made with natural rubber latex
- Smart Dispense Technology
- Tested for protection against chemotherapy drugs
- · ESD Tested and approved

## **BENEFITS**

- · Strong and more durable
- Easy donning, improved fit and comfort
- Eliminates exposure to latex
- Single glove delivery with fewer wrinkles
- Appropriate for handling cytotoxic chemicals
- · Safe for use with sensitive electronics

## PRODUCT CHARACTERISTICS

- Synthetic Nitrile Rubber (NBR)
- Tested for Use with Chemotherapy Drugs
- Textured Fingers & Palm
- Powder-Free
- Ambidextrous
- **Beaded Cuff**
- Single Use Only



orth a Tra
SPROVED * TESTED
STED.
日
FSD.

Order#	Size	Gloves/Box*	Boxes/Case
NIPFT101	XS	100	10
NIPFT102	S	100	10
NIPFT103	М	100	10
NIPFT104	L	100	10
NIPFT105	XL	100	10

\*By weight

Gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. Review material safety data sheets for the drug being used to determine the required level of protection. Test results available upon request. Chemotherapy drugs tested.



## HOW WE GLOVE THE PLANET WITH SEMPERCARE® NITRILE

✓ Inner boxes made with 80% recycled material Outer cartons made with 90% recycled material



www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762 Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558



Glove the Planet is Sempermed USA Inc.'s commitment to reducing our environmental impact. By thinking ahead and consuming less we're able to do more to minimize our global footprint with every step we take. Our ethos of sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations. To learn more, visit www.GlovethePlanet.com.

