



Gel Tubes

Silicone gel tubes are ideal to cushion fingers or toes.





118 Log Canoe Circle
Stevensville, MD 21666 U.S.A.

TO ORDER

Tubes stretch to fit most fingers or toes.

One size 1/2" x 1 3/4" (1.5cm x 4.5cm)

P1111 - Package of 5

Order Toll Free 888.378.7763 (U.S only)

Local/Int. 410.604.6393

24 hour Fax 410.604.6398

Mon. - Thur. 8:30 - 5:30 pm, Fri. 8:30 - 5:00 pm EST

CARE INSTRUCTIONS

Wash in cold water with mild soap or detergent. Tubes may darken with wear. Change in color will not affect the product's performance.

Material: Silicone gel. Contains no latex.

www.3pointproducts.com

Authorized European Representative
Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands




©2009, 3-Point Products, Inc.

Made in Taiwan


Cushioning Oval-8® Finger Splints with Gel Tubes

Gel Tubes from 3-Point Products are ideal to cushion finger splints including 3-Point Products' Oval-8® Finger Splints. The instructions below demonstrate cutting a Gel Tube to cushion Oval-8 splints with little or no change in splint sizing.

Cutting instructions to create an "apron".



Gel Tube used to cushion the Oval-8 band. Shown for a Mallet Finger application.



Gel Tube to cushion one or both ends of an Oval-8 splint. Shown for a Trigger Finger application.

