



Stainless Steel Emesis Basin 6 12oz ⁴²³⁰

12/BX

Product Features	Made of the high grade stainless steel (18/8,304)	
	Crafted without seams which prevents bacteria growth and allows for easy cleaning	
	Autoclavable	

Attributes	Brand Name	Tech-Med
	Size	6"
	Material	18/8 Stainless Steel
	Sterile	No
	Weight	3.04238 lbs
Packaging	Packaging	12/BX
	Case Packaging	12/BX
	Total Case Quantity	12
	Item Packaging Barcode	665973021737
	Inner Packaging Barcode	(01)10665973021734
	Case Packaging Barcode	(01)20665973021731

Phone: 1 (631) 656-3800 **Toll Free:** 1 (800) 243-0741 **Fax:** 1 (631) 656-3810

