

5-Leg Pneumatic Stool

MODEL 2155

- Premium, 23" die cast, aluminum base with satin finish top
- 5-leg design for greater stability
- Full size diameter top and base for comfort and stability
- 4" thick seat for extra comfort
- Black hooded, dual wheel, nylon casters
- Pneumatic height adjustment
- All seats have 3/4" thick solid plywood base
- Lever controls height
- Knock-down design saves on shipping
- Easy, no-tools required assembly
- Premium, stain-resistant, woven, knit-backed vinyl upholstery
- Also available with backrest
- 250 lbs. load capacity under normal use



Seat Diameter	16"
Base Diameter	23"
Height Range	19"–24"



California residents please be advised of the following, pursuant to Proposition 65:

warning: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Options



Upholstery Colors

