



Sterile Wound Closure Strip $1/4 \times 3$

5151

5/PK, 50 PK/BX, 4 BX/CS

Product Features

Adhesive backed, thin non-woven, reinforced fibers that are arranged in a bilateral pattern for high tensile strength

Provides support following early suture/staple removal

Keeps wound closed without significant tension, edema, stretch or movement

Translucent and porous

Attributes

Brand Name	Dukal
Size	1/4 x 3
HCPC Code	A4450
Sterile	Yes
Weight	1.3 lbs

Packaging

Packaging 5/PK, 50 PK/BX, 4 BX/CS

Item Packaging 3/PK **Inner Packaging** 50 PK/BX **Case Packaging** 4 BX/CS **Total Case Quantity** 600

Item Packaging

(01)10665973008841

Inner Packaging

Barcode

Barcode

(01)20665973008848

Case Packaging

Barcode

(01)30665973008845

2 Fleetwood Ct. Ronkonkoma, NY 11779 Phone: 1 (631) 656-3800 **Toll Free:** 1 (800) 243-0741 **Fax:** 1 (631) 656-3810

