

# XL CONVALESCENT RECLINER

5291 (WITH DUAL SIDE-TABLES)

5 SERIES



## Superior Comfort and Support...

- Built-in headrest
- Removable side panels for ease of cleaning
- 3 Locking positions: upright, TV or recline
- Dual side tables
- Padded armrests with scuff guards
- Recessed handgrips provide additional assistance for ingress and egress


ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



## SPECIFICATIONS

DROP-ARM CONVALESCENT RECLINER	5291
DESCRIPTIONS	DUAL SIDE TABLES INCLUDED
Seat Height	22.5" (54.2 cm)
Seat Width	25" (63.5 cm)
Seat Depth	20" (50.8 cm)
Overall Height	49" (124.5 cm)
Overall Width (side tables down)	36.25" (90.1 cm)
Range of Seat Angle	8° - 32°
Range of Back Angle	62° - 38°
Back Height above Seat	31.5" (80 cm)
Overall Length (Back Upright)	44" (111.8 cm)
Overall Length (Back Reclined)	62" (157.5 cm)
Weight	112 lbs (50.8 kg)
Weight Capacity	450 lbs (204.1 kg)

## SUGGESTED OPTIONS & ACCESSORIES



**IV POLE**  
Left, Right, or  
Dual

**OVER 50  
DESIGNER  
COLORS  
AVAILABLE**

*Additional options and accessories available.*

*Please visit our website to view a full selection of options or contact your local representative.*

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

WMF5291 REV B