

Reorder No. 5627

Cardboard Head Immobilizer



- 1 Made of sturdy moisture-resistant cardboard.
- 2 Includes comfortable foam head pad.
- 3 Quickly attaches to board with EZ peel tabs, no straps or clips needed.
- 4 Radiolucent.
- 5 Patient's ears are easily accessible.
- 6 Universal and adjustable – will fit virtually all backboards and patients.
- 7 Individually sealed in plastic for protection from the elements until ready to use.
- 8 Compact size—stores flat.
- 9 Disposable. Intended for single use only by authorized EMS personnel.
- 10 Writable surface.



 **Manufactured for:**
Dynarex Corporation
10 Glenshaw Street
Orangeburg, NY 10962
USA • www.dynarex.com

Made in China



SYMBOL GLOSSARY
For an explanation of symbols
used in Dynarex packaging, visit
dynarex.com/symbols.php

**For More Emergency
Medical Service Products,
go to www.dynarex.com.**

© 2018 Dynarex Corporation.

R181112