Basic Quick-Release Head Immobilizer dynarex







- Compatible with all C-collars and backboards
- Radiolucent
- Easy access to the ears
- Head strap sticks to board, not to skin/hair (foam covers forehead area)
- Non-claustrophobic design
- · Closed-cell foam makes it water resistant and buoyant for water rescue
- · Head and chin straps included with each unit
- Can be used with motorcycle and football helmets

How to use:







Manufactured for: **Dynarex Corporation** 10 Glenshaw Street Orangeburg, NY 10962 USA • www.dynarex.com Made in China

Cancer and Reproductive Harm www.P65Warnings.ca.gov

For More Emergency Medical Service Products, go to www.dynarex.com

© 2018 Dynarex Corporation R180913