

Ocean Commotion, 2-Bin, Phlebotomy Cart

MODEL 67236

- Clinton Kids™ Ocean Commotion graphics embedded in laminate will not tear or peel
- Lockable door on soft close hinges
- Two, 3.5" x 11" x 16.5" removable plastic storage bins on easy slide tracks (9 cm x 28 cm x 42 cm)
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Easy grip, ClintonClean™ solid plastic handle
- Pneumatic height adjustment
- Black, nylon 5-leg base
- 2.5" dual wheel casters
- Optional glove box holder available. Add product GW-2000 Single White Steel Glove Box Holder



NOTE: Items pictured on top not included

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

MODEL 67236 SPECIFICATIONS

Top	20" W x 20" D
Base Diameter	24"
Height Range	33" — 43"
Cabinet	13.5" W x 18" D x 10.5"



Prop 65 Warning



WARNING: This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: www.P65Warning.ca.gov.

7

OPTIONS



GW-2000 Single White Steel Glove Box Holder