



Completely Coordinated.
Completely Cost Effective.
Completely Easy.

ZOO BUS COMPLETE



Fun Series, Zoo Bus, Treatment Table



FRONT & DOOR SIDE



BACK & OPPOSITE SIDE

MODEL 7020

- Clinton original illustrations
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on two long side and two short ends, manufactured as pictured
- Optional locking doors available
- All colors are as shown
- 3 Storage compartments with one adjustable shelf each
- Soft-close, door hinges
- Chip resistant PVC edges
- Seamless rounded corner top
- Premium, stain resistant, woven knit back upholstery (black)
- Adjustable feet
- Doors are only on one side
- Paper dispenser included
- Special order reverse side doors are available
- 350 lbs. load capacity under normal use

7020 SPECIFICATIONS

72"L x 27"W x 34"H

California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Options




066 Door & Inside Latch Combo




055 Door Locks



040 Paper Cutter



20 Small Pillow



60 Large Pillow

Upholstery



Black

Schoolhouse Cabinets

MODEL 6120-BW

"Schoolhouse" Base Cabinets

4 doors with adjustable shelves and 3 drawers.
Shipped in one unit.

"Schoolhouse" Wall Cabinets

4 doors with adjustable shelves. Shipped in one unit.

Base & wall cabinets are only sold in sets

- Clinton original illustrations
- Pre-assembled units save time
- Easy-clean all laminate construction
- All graphic are permanently embedded in the high pressure laminate and will not peel off
- Tough, chip-resistant PVC edges
- Assembled with steel fasteners
- 4 adjustable leveling feet on each unit
- 3-way adjustable, soft-close 110° hinges
- Flat panel drawer fronts with metal sides that hold up to 75 lbs.
- Recessed toe kick
- Easy hook and rail installation for upper cabinets
- Slate Gray edgebanding on all edges
- Customizable with 4 options
- Easily moves for relocation
- Drawers can be added to any base unit
- Doors and drawer can be moved to accommodate sink location
- Shown with optional sink and faucet



SPECIFICATIONS

BASE: 60"L x 18"D x 35"H


WALL: 60"L x 12"D x 24"H

California residents please be advised of the following, pursuant to Proposition 65:


! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Options



066 Door & Inside Latch Combo



055 Door Locks




044 Add a Cabinet Drawer




22 Sink & Faucet

Optional Items in Red



1 Square = 1 foot

NOTE: These are suggested floor plans. Please check local building and fire codes before beginning any construction or renovations.



Paint option with chair rail



Walls

2022-40 Banana Yellow

795 Highland Blue

Chair Rail/Trim

797 Athens Blue

Flooring

Medintone® H5301 Gray Light

Medintone® H5302 Deep Gray

OPTIONAL ACCESSORIES *(Shown in photo on front)*



- 01-CC Koala Pyramid Wall Sticker (34" x 73")
- 10-K Koala Stool
- C-40 Side Chairs in Royal Blue upholstery (3RB)
- T-70 Clear Labeled Sundry Jars (5 in set)
- GS-3022 Double Slimline Stainless Steel Glove Box Holder



800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

www.clinton-ind.com