





Re-Order#	Description	Size	Unit
703060	Leg Bolster	4"D x 6.5"H	Pair
703061	Leg Bolster	4"D x 6.5"H	16 Pairs/CS

Product Description

Skil-Care Leg Bolsters is designed to guard against skin tears due to sharp and/or protruding edges around footrest/legrest hanger brackets.

Application

Wrap the Skil-Care Leg Bolster around the footrest/legrest hanger bracket and secure with hook & loop fastener strap.

Easy Care

The LSII cover can be wiped clean using a damp cloth, spray-type cleanser or mild disinfectant. Do not use bleach.