

Fun Series, Zoo Bus, Pediatric Scale Table



MODEL 7822

TABLE FEATURES

- Clinton original illustrations
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Built-in concealed paper dispenser and cutter
- Adjustable built-in levelers provide a stable surface
- Easy-clean all laminate base
- Optional storage area locks available
- Ample 2-door storage
- Soft-close, door hinges

SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.



MODEL 7822 SPECIFICATIONS

Top:	43"L x 24"W
Base:	47"L x 21"W x 36"H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.




Upholstery



Gunmetal

Features



Infantometer

Built-in infantometer with folding ends for measuring height



Padded Top

Easy-clean, padded, upholstered top for comfort



Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



Scale Head

Easy-read, digital display