


## Imaging Table with Fowler Back and Drop Window



### MODEL 80072

- Two-piece top with pneumatic adjustable backrest and window drop
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- 34" wide model available
- Model 80072: 500 lbs load capacity under normal use
- Model 80072-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



*Trendelenburg*



*Reverse Trendelenburg*



*Optional Side Rails & Casters*

### 80072 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 <sup>1</sup> / <sub>2</sub> " – 37 <sup>1</sup> / <sub>2</sub> "
backrest	38"	27"	
window	10"	10"	

**80072-X** Extra-Wide 34" width model



### Upholstery Colors



## OPTIONS

### Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm

### Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098

### Power Backrest

- Smooth, power control raises to 80° and lowers to flat

011

### Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

### Paper Cutter

- Clear plastic with elastic strap and buckle

040

