seca aura 807

The first picture frame for the floor.



- Four-cell weighing technology. For extra flat construction.
- Step-off function. For weighing without switching on scales.
- Safety glass. For extra big impact and shock resistance.
- 1.1 inch LCD figures. For easy reading.
- More than 165 years of experience. For maximum precision.
- From the global market leader for medical weighing and measuring.



seca aura 807:

Projecting a beautiful image.

ake a look at the design: the display is not just a display, but a work of art in a frame of brushed stainless steel. This extra is only one of many details which speak for the class of the seca aura 807.



Projecting a high-quality image.


It not only looks a bit like a picture frame, it also looks like it should be in one. This attractive scale has an extremely flat weighing platform of single-pane safety glass, an extra flat base of high-quality thermoplastic and four elegant load cells of brushed stainless steel. All this makes it exactly the right scale for people who appreciate beauty.

Projecting a timeless image.

The elegance of the design does full honour to its name "aura", a fact underscored by the 1.1 inch LCD figures in a frame of brushed stainless steel. Be reassured - the glass plate is out of harm's way, as it is somewhat smaller than the frame below. The best view isn't from above but from the side. Check out the extremely flat body with its foot-friendly, rounded edges.

Giving everything, plus a little more.

There's also something extra that's only available with personal scales from seca: the Step-off function. You can use it without switching on the scale. Simply stand on the platform, wait for the "Step-off" instruction, and then Step off to the side and read off the weight in the display. You see, a personal scale from seca isn't just good, but in many aspects simply better.



Display frame made of brushed stainless steel.



seca aura 807

Technical data

- Capacity: 330 lbs / 150 kg
- Division: 0.2 lbs / 100 g
- Dimensions scale (WxHxD): 12.2 x 1.5 x 12.2 inch
- Weight: 4.2 lbs
- Functions: Step-off, Auto-HOLD, automatic switch-off
- Parts to be delivered: 4 AAA batteries, 4 extension feet for carpets