


Open Base, Extra Wide, Bariatric, Straight Top Power Table

MODEL 84100

- Extra wide, bariatric model with choice of two widths
- Features easy-clean, one-piece top
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 600 lbs. (272 kg) load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:



POWER 600 BARIATRIC



84100-34 SPECIFICATIONS

Overall

72" L x 34" W x 18-35" H

84100-40 Extra-Wide 40" width model





OPTIONS

Safety Straps Set

- 2"x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 safety straps per set (92" x 2")

080

Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082

Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084

Locking Casters

• Heavy duty, rubber wheel, easy-roll casters

087

030 Paper Dispenser

• Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

Paper Cutter

• Clear plastic with elastic straps and buckle

040

Paper Dispenser/Cutter Combo

• Paper dispenser and paper cutter ordered as one option

034





