

E-Series, Power Table with Adjustable Backrest and Drop Section

MODEL 86300

- Features an adjustable, pneumatic backrest and friction lock drop section which are accessible from either side
- Solid, standard drop section (may be ordered with a nose hole as an option) (Please specify when ordering)
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 1.5" (4 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 450 lbs. (204 kg) load capacity under normal use

E-Series
POWER TABLES



SPECIFICATIONS

Overall 76" L x 27" W x 18.5 –35.5" H

Backrest 38" L x 27" W

Drop Section 13" L x 27" W


California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Options



009 Nose Hole — Nose hole in drop section of 86220 or 86300

030 Paper Dispenser — Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

040 Paper Cutter — Clear plastic with elastic strap and buckle

034 Paper Dispenser/Holder Combo — Paper dispenser and paper holder ordered as one option

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black