

MONARK 881E

MAIN FEATURES

- Arm and leg ergometer
- Individually adjustable crankarms
- Possible to calibrate

THE BASICS

- Good cycle feeling
- User friendly display
- Can be used from wheelchair
- No power adaptor needed
- Pedals and handles included

MONARK REHAB TABLE

- Adjustable hight
- Adapted for 881E
- Lockable
- Wheelchair size

Medical Device class lla



Specifications

MONARK 881E

Certifications	EU directive
Cei tilications	93/42/EU class lla
Frame	
Туре	Arm ergometer
Handlebar	n/a
Stem	n/a
Seatpost	n/a
Saddle	n/a
ВВ	n/a
Crank	Steel, 152 mm
Pedals	9/16", handles or 9/16" pedals with foot straps
RESISTANCE	
System	Spring
Flywheel	11 kg
DISPLAY	
RPM	Yes
HR	n/a
Time	Yes
Speed	n/a
Distance	n/a
Watt	n/a
Connectivity	
RS232	n/a
Usb	n/a
WEIGHT & MEASUREMENTS	
Max user weight	n/a
Length (mm)	540 mm
Width (mm)	470 mm
Height (mm)	550 mm
Weight (kg)	22 kg
POWER REQUIREMENTS	Batteries in the display



WHEELCHAIR

Adapted for use from a wheelchair

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 15 36 01 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden

Upplaga 2. 2018-06-19

