

## Family Practice Ready Room



*5* Year  
Limited  
Warranty

CLINTON  
**Ready Room**™

*an all-inclusive Clinton exclusive*

A complete exam room, well equipped and ready for patients. Clinton's 8890 Family Practice Ready Room package offers big advantages in function and added accessories at a much smaller price point than other manufacturers.

*(no model substitutions, your choice of Classic Clinton laminate and upholstery colors)*

# MODEL 8890-RR

## PACKAGE INCLUDES:

- 8890 in your choice of laminate and upholstery colors
- 8248 Wall Cabinets in your choice of laminate colors
- 8048-022 Base Cabinets in your choice of laminate colors with stainless steel sink and chrome faucet set included (choose right or left location or no sink. If no location is selected, sink will be placed as shown.)
- 21335 in your choice of upholstery colors
- GS-3000 Single, Stainless Steel,\* Glove Box Holder
- TR-20S 20-Quart, Medium Round, Stainless Steel,\* Waste Receptacle
- C-40B Premium, Padded, Side Chair in your choice of upholstery colors

\* Not MRI Compatible

## 8890 Family Practice Table with Step Stool

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 28.5" width for patient comfort (72 cm)
- Handy, pull-out, welded steel, step stool
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Easy care all laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 drawers with all steel slides
- Built with plywood and laminate panels that are California CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181 kg)



### SPECIFICATIONS

58"/72" extended L x 28.5" W x 31" H

### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### Laminate Colors



Dark Cherry



Maple



Gray



Slate Gray

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

## 8248 Wall Cabinet

- Wall cabinet with 2 doors and 1 adjustable shelf
- Easy-clean, all laminate surface
- Tough, chip resistant edges
- Unit is pre-assembled to save installation time
- Cabinet easily mounts with hook and rail mounting system
- Satin finish pulls
- Choice of solid color or woodgrain laminates
- Fully adjustable, concealed, soft-close hinges

*fashionfinish*  
CASEWORK BY CLINTON



## 8048 Base Cabinet

- Base cabinet with 2 drawers and 4 doors with 2 adjustable shelves
- Stainless steel sink and chrome faucet (left or right placement)
- Easy-clean, all laminate surfaces
- Tough, chip resistant edges
- Unit is pre-assembled to save installation time
- Fully adjustable, soft-close hinges
- 3.25" backsplash included with top
- Satin chrome pulls
- Choice of solid color or woodgrain laminates



### 8248 SPECIFICATIONS

48" W x 12" D x 24" H

### 8048 SPECIFICATIONS

48" W x 18" D x 35" H

### Prop 65 Warning

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Finishes



Maple



Dark Cherry



Gray



Slate Gray

## 2135 Standard Series, 5-Leg, Pneumatic Stool

- Thick padded seat for comfort
- 24" Black nylon, 5-leg base
- Nylon dual wheel casters
- Easy, pneumatic height adjustment
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- Maximum load capacity 250 lbs. under normal use



### SPECIFICATIONS

Seat Diameter 16"

Base Diameter 24"

Height Range 19"-24"

### Prop 65 Warning

**!** **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

### GS-3000 Single, Stainless Steel,\* Glove Box Holder

- Stainless steel construction
- Open design allows for easy glove identification
- All units mount horizontal or vertical
- Mounting hardware included
- Universal size fits most glove boxes

\*Clinton stainless steel items are not M.R.I. compatible.



Shown in the horizontal position

#### SPECIFICATIONS

10" L x 3.75" D x 5.5" H

#### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### TR-20W Medium, Round, White, Waste Receptacle

- All steel outer shell and lid
- White enamel finish
- Step-on foot pedal operated lid
- One-piece, plastic base keeps metal from contact with floor surface
- Chrome Handle make moving and cleaning easy
- 2 warning labels included
- 20 L capacity (21.13 quartss)



#### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

#### SPECIFICATIONS

10" D X 24"H

### C-40B Premium, Padded, Black Frame Side Chair Without Arms

- Black powder coated frame finish
- Heavy duty 1-1/2" diameter legs
- Solid plywood seat and back
- Wall saver guard
- 2" thick, firm foam padded seat and back
- Contoured seat and back
- Adjustable leg levellers
- Premium, knit-backed upholstery
- 300 lbs. load capacity under normal use



#### Prop 65 Warning

**WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: [www.P65Warning.ca.gov](http://www.P65Warning.ca.gov).

7

#### SPECIFICATIONS

Seat 18" W x 18" D

Seat Height 18"

#### Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black