

Fashion Finish Family Practice Ready Room



CLINTON
Ready Room[™]

an all-inclusive Clinton exclusive

Get Fashion Finish casework packaged with 5 additional, essential items for the ultimate in modern exam room function and decor. And best of all, it's easy to order, durable enough for years of service.

(no model substitutions, your choice of Clinton Fashion Finishes, countertop colors and upholstery colors)

MODEL 8890-RRFF

PACKAGE INCLUDES:

- 8890 in your choice of laminate and upholstery colors
- 8348 Wall Cabinets in your choice of Fashion Finishes
- 8648-022 Base Cabinets in your choice of Fashion Finishes with stainless steel sink and chrome faucet set included (choose right or left location or no sink. If no location is selected, sink will be placed as shown.)
- 48P postform countertop in your choice of 3 colors
- 21335 in your choice of upholstery colors
- GS-3000 Single, Stainless Steel,* Glove Box Holder
- C-40B Premium, Padded, Side Chair in your choice of upholstery colors

* Not MRI Compatible

8890 Family Practice Table with Step Stool

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 28.5" width for patient comfort (72 cm)
- Handy, pull-out, welded steel, step stool
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Easy care, Gray, laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 drawers with all steel slides
- Built with plywood and laminate panels that are California CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181 kg)



SPECIFICATIONS

58"/72" extended L x 28.5" W x 31" H

Prop 65 Warning

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Laminate Color



Gray

Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black

8348 Wall Cabinet

- Wall cabinet with 3 doors and 1 adjustable shelf
- Durable, easy-clean, thermofoil surfaces
- No front facing seams
- Unit is pre-assembled to save installation time
- Cabinet easily mounts with hook and rail mounting system
- Fully adjustable, concealed, soft-close hinges
- Smart looking, easy-clean, satin chrome pulls
- Choose from 5 contemporary finishes

fashionfinish
CASEWORK BY CLINTON



8648-022 Base Cabinet

- Base cabinet with 3 doors, 2 drawers and 2 adjustable shelves
- Stainless steel sink and chrome faucet set
(Position left or right)
- Easy-clean, thermofoil surfaces
- No front facing seams
- Units are pre-assembled to save installation time
- Fully adjustable, soft-close hinges
- Smooth glide, drawers with one-piece metal sides hold up to 75 lbs.
- Smart looking, easy-clean, satin chrome pulls
- Choose from 5 contemporary finishes



8348 SPECIFICATIONS

48" W x 12" D x 24" H

8648 SPECIFICATIONS

48" W x 18" D x 35" H

Prop 65 Warning

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Finishes



Chestnut Hill



Twilight



Arctic White



Metropolis Gray

48P Postform laminate countertop

Postform laminate countertops offer a distinctive one-piece countertop that eliminates the seam between countertop surface, backsplash and front edge.

A stainless steel sink and chrome, wing lever faucet is standard with all postform countertops. The faucet is 14.5" high with a 3" radius.

Sink may be positioned on the left or right, or center on a 72" cabinet. Please specify when ordering base cabinet and top.

Model includes sink and faucet. All models with postform countertops without sinks are special order.



SPECIFICATIONS

48P—48"L x 19"D X .75"H


Backsplash —3"

Prop 65 Warning

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Postform Colors



White Carrera



Meteorite



Black Alicante

2135 Stand Series, 5-Leg, Pneumatic Stool

- Thick padded seat for comfort
- 24" Black nylon, 5-leg base
- Nylon dual wheel casters
- Easy, pneumatic height adjustment
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- Maximum load capacity 250 lbs. under normal use

SPECIFICATIONS

Seat Diameter	16"
Base Diameter	24"
Height Range	18"-23"



Prop 65 Warning

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

GS-3000 Single, Stainless Steel,* Glove Box Holder

- Stainless steel construction
- Open design allows for easy glove identification
- All units mount horizontal or vertical
- Mounting hardware included
- Universal size fits most glove boxes

*Clinton stainless steel items are not M.R.I. compatible.

SPECIFICATIONS

10" L x 3.75" D x 5.5" H



Shown in the horizontal position

Prop 65 Warning

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

C-40B Premium, Padded, Black Frame Side Chair Without Arms

- Black powder coated frame finish
- Heavy duty 1-1/2" diameter legs
- Solid plywood seat and back
- Wall saver guard
- 2" thick, firm foam padded seat and back
- Contoured seat and back
- Adjustable leg levers
- Premium, knit-backed upholstery
- 300 lbs. load capacity under normal use

SPECIFICATIONS

Seat	18" W x 18" D
Seat Height	18"



Prop 65 Warning

! **WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: www.P65Warning.ca.gov.

7

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black