

## Style Line, Cabinet Style, Treatment Table with Sliding Doors



### MODEL 9044

- Features an adjustable backrest, 2 storage compartments, 2 sliding doors and an adjustable shelf in each compartment
- Easy-clean laminate surfaces
- Ships assembled
- Seamless rounded corner top
- Heavyweight, knit-backed upholstery
- 2" firm foam padding (5 cm)
- 425 lbs. load capacity under normal use (193 kg)

**Note: If selecting 055 door locks, 2 door locks MUST be ordered to secure the storage compartment.**

### SPECIFICATIONS

9044-27	72" L x 27" W x 31" H
9044-30	72" L x 30" W x 31" H

**California residents please be advised of the following, pursuant to Proposition 65:**

**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



### Laminate Colors



### Upholstery Colors

