

Style Line, Cabinet Style, Treatment Table with 4 Doors



MODEL 9074

- Features 4 doors, adjustable backrest and 2 storage compartments each with an adjustable shelf
- Easy-clean laminate surfaces
- Adjustable doors with concealed soft-close hinges
- Ships assembled
- Seamless rounded corner top
- Heavyweight, knit-backed upholstery
- 2" firm foam padding
- 425 lbs. load capacity under normal use


SPECIFICATIONS

9074-27	72" L x 27" W x 31" H
9074-30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Options



Laminate Colors



Upholstery Colors

