

Fetal Monitor Cart

Model 9202-C



**shown with optional
monitor mount*


Features:

- Constructed of medical grade Melamine material
- Pull Out Tray, 2 drawers and bottom compartment
- Cord Cut-Out to allow for easy use as a monitor stand
- Extender Shelf: 14"
- 5" locking casters
- 1 Year Limited Warranty

Overall Dimensions

42"H x 22"W x 19"D

Finish Options:



Prestige Maple Candlelight Classic Cherry Biltmore Cherry Truffle

Medical grade case carts constructed of aesthetically pleasing melamine material to create a warm medical environment. Designed to meet all of your practical needs with their various compartments and drawers