

NYOrtho

NYOrtho ZERO-G Boot™

MODEL 9518

The NYOrtho ZERO-G Boot™ helps eliminate pressure and friction to the skin and floats the heel over an air cavity providing air circulation to aid the healing of existing ulcers.

RECOMMENDED USE: Individuals unable to avoid long periods of uninterrupted pressure over bony prominences are at increased risk for the development of necrosis and ulceration.



ANTIMICROBIAL

No odor, no bacteria: thanks to the built-in hygiene function.



SUPER-SOFT VELVET

Protects thin fragile skin from tears caused by trauma, reducing friction and shear.



MACHINE WASHABLE

View care label for laundering instructions.



ADAPTABLE STRAP

Permits this device to be secured at two alternate positions. Lower position controls dorsiflexion.

NYOrtho ZERO-G Boot™

Illustrated with cutout showing heel floating over air cavity, foam layers and gel-pack.



PADDED HOOK AND LOOP

Easy-on, easy-off padded closure allows easy application and removal with a comfortable fit.

FOAM CUSHIONED SIDES

Lateral and medial malleoli – bony prominence on both sides of the ankle, are softly cushioned.

ZERO PRESSURE

Extended cushioned design floats the heel over an air cavity for ZERO pressure and provides air circulation to aid the healing of existing ulcers.

A REMOVABLE GEL-PACK

May be utilized to further lower pressure beneath the ankle and keeps feet cool. (FA-GEL-4X8)

HOW TO MEASURE

To choose the correct size, measure the circumference of the lower part of the calf muscle.

Model	Size	Circumference	Fits
9518-S	Small	5-10 in.	Petite
9518-M	Medium	11-20 in.	Adult
9518-L	Large	21-25 in.	Bariatric

Independent laboratory tests via C-PicoPress Interface Pressure Monitor, with ambient temperature consistent with ISO 554-1976E, prove pressure is effectively removed from the heel. EC Service Project#1385

