



Informed  
conversations  
at the  
point of care

## HbA1c results in 5 minutes

- Improve patients' in-office experience
- Increase efficiency
- Help drive better health outcomes for diabetes patients

**A1CNow<sup>+</sup>**



### A1CNow<sup>+</sup> Test System Key Features

- Lab-quality accuracy
- Portable
- Easy to use
- Fingerstick (5  $\mu$ L) blood sample
- NGSP-certified
- IFCC-traceable
- CLIA-waived

# A1CNow<sup>+</sup> Controls™ and Connectivity



## Ensure Accurate A1C Results



**A1CNow<sup>+</sup> Controls™**  
REF 737

- Ready-to-use control materials
- Two-level control kit
- Dropper vials



**Room Temperature Storage**

7 Days  
once thawed and opened  
18° to 25°C (64° to 77° F)

**Refrigerated Storage**

8 Months  
once thawed and opened  
2° to 8°C (36° to 46° F)


**Long-Term Frozen Storage**

3 Years  
or until expiration date  
-30° to -10°C (-22° to 14° F)

## PTS Connect Blu-Dock™

*(Sold separately)*

REF 4715



- Wirelessly transfers results to third-party data management systems<sup>1</sup>
- Bluetooth enabled
- Seamlessly connects with a PTS Connect Printer

## PTS Connect Printer

*(Sold separately)*

REF 2750



- Outputs results on paper or self-adhesive labels
- Can connect to A1CNow<sup>+</sup> through Blu-Dock

Available in a 10-count or 20-count test kit



**A1CNow<sup>+</sup> system 10-count test kit**  
REF 3024



**A1CNow<sup>+</sup> system 20-count test kit**  
REF 3021

<sup>1</sup> May require additional development.

**For more information or to place an order:**



Website: [ptsdiagnostics.com](http://ptsdiagnostics.com)  
Office: +1 (317) 870-5610 FAX: +1 (317) 870-5608  
Toll-Free in the US: 1 (877) 870-5610