



EssentialCARE Anti-Embolism Stockings provide classic comfort and features with an integrated color-coding system for easy size identification.

Alignment Indicator

Two-Way Stretch Fabric

Color-Coded Heel Positioner

High visibility size identification. Assists clinician with proper positioning of stocking.

Top View Inspection Port

## EssentialCARE® Features and Benefits

### Patient Comfort and Safety

**Top-View Inspection Port** Provides better patient comfort when ambulating.

**Rapid Moisture Transfer** Moves fluids and body oils from skin to surface of the fabric quickly to speed drying; launders better.

**Two-Way Stretch Fabric** Allows stocking to conform to shape of limb.

**Premium Elastic Top** Keeps stocking in place without causing tourniquet effect.

**Outward Facing Seams** Reduces skin irritation and/or depressions as limb size changes.

**Latex-Free** Safe for latex-sensitive skin.

### Clinician Focused

**Color-Coded Heel Positioner** Assists proper heel and ankle placement; easy size identification.

**Color-Coded Top-View Inspection Port** Conveniently access toes and top of foot and identify size.

**Alignment Indicator** Visible guide to prevent twisting during application.

**20-30mmHg Graduated Compression** Helps improve blood flow and aids in clot prevention.

### Available Sizes and Styles


**32 Sizes To Choose From** Knee and thigh lengths.

For more information, visit [www.DVTSTAR.com](http://www.DVTSTAR.com)



## EssentialCARE® Anti-Embolism Stockings

Short	Regular	Long	Style	Size	Color	Qty.
951-01	953-01	958-01	Knee	Small	Yellow	12 PR
951-02	953-02	958-02	Knee	Medium	Orange	12 PR
951-03	953-03	958-03	Knee	Large	Blue	12 PR
951-04	953-04	958-04	Knee	X-Large	Green	12 PR
951-05	953-05	958-05	Knee	2X-Large	Grey	12 PR
—	953-06	—	Knee	3X-Large	Purple	12 PR
963-01	968-01	973-01	Thigh	Small	Yellow	12 PR
963-02	968-02	973-02	Thigh	Medium	Orange	12 PR
963-03	968-03	973-03	Thigh	Large	Blue	12 PR
963-04	968-04	973-04	Thigh	X-Large	Green	12 PR
963-05	968-05	973-05	Thigh	2X-Large	Grey	12 PR
—	968-06	—	Thigh	3X-Large	Purple	12 PR



Encompass Group, LLC  
615 Macon Street  
McDonough, Georgia 30253  
800-284-4540  
[www.encompassgroup.net](http://www.encompassgroup.net)

 **encompass®**  
Better care starts with safety and comfort.