

ENCORE[®] Latex Acclaim

S U R G I C A L G L O V E S

Latex, powder-free surgical gloves, offering wearers a comfortable, natural fit and feel

- **Superb dexterity and tactile sensitivity:** These latex surgical gloves' optimal thickness, made possible due to their soft, yet robust formulation, offers wearers the dexterity and tactile sensitivity needed for safe and efficient performance
- **Elevated ergonomic comfort:** Together, the unique anatomic formers and soft latex formulation of these powder-free surgical gloves have resulted in certification for ergonomic comfort by United States Ergonomics
- **Excellent grip:** ENCORE[®] Latex Acclaim surgical gloves feature a smooth external surface, with a micro-textured finish, enabling the safe and secure grip of medical instruments especially when conditions are wet
- **Enhanced donning aid:** Ansell's proprietary donning aid technology ensures effortless donning, be it the first don or during intraoperative changes
- **Tested for chemo handling:** These latex surgical gloves have been tested in accordance with ASTM D6978 and may be safely used to handle a range of chemo drugs*



CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- Tested for use with chemotherapy drugs
- Extra comfort during long, sustained procedures
- Damp donning
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES



ENCORE[®] The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

ENCORE[®] Latex Acclaim

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D3577, ASTM D6978, EN 16523-1, EN 374-2 and -4, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	5795000	5795001	5795002	5795003	5795004	5795005	5795006	5795007

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com