Baxter

I·I[•] Hillrom_•





SURETEMP PLUS 692 SETTING THE STANDARD OF CARE

Advance your temperature routine with the **SureTemp** Plus 692 Thermometer for hospitals, featuring technology that sets the standard for speed and accuracy.

1 EASY

Intuitive user interface, easy-to-read LCD screen and flashing icons, helps select and identify temperature measurements and modes.

2 SAFE

Designed with interchangeable oral/axillary (blue) and rectal (red) probe wells to address patient safety and help reduce the risk of cross-contamination.

3 FAST AND ACCURATE

With proven accuracy, the **SureTemp** 692 thermometer captures quick readings-approximately 4-6 second oral, 10-13 second pediatric axillary and rectal readings and 12-15 second adult axillary readings.

4 SECURE

To top it all off, our security features make sure the unit stays where it is supposed to be with a location ID on the display and other innovative features help prevent theft. Interchangeable probe with push-button probe cover ejection

Removable probe well for easy cleaning

C/F button for toggling between Fahrenheit and Celsius

Prominent recall button for display of last temperature taken

Ergonomic design with rubberized grips

Welch Allyn SureTemp Plus 692

2 WES1

0

User-selectable icon and/or words for identification of thermometer mode Convenient storage for 25 probe covers

Large, easy-to-read, backlit LCD

Mode button for selecting oral, pediatric axillary or adult axillary

Pulse timer button/ monitor mode

NEW feature makes it easier to switch to Monitor Mode for measuring a patient's thermal steady state (approx. 3 to 5 minutes)

Interchangeable rectal probe and well reduce risk of cross-contamination

annun nu

STANDARD WALL MOUNT

The **SureTemp** Plus 692 with standard wall mount has storage for an additional 50 probe covers.



Welch Allyn **SureTemp** Plus 692 Electronic Thermometer

SETTING THE STANDARD OF CARE EVEN HIGHER

TECHNICAL SPECIFICATIONS

Calibration Accuracy	±0.2 °F (±0.1 °C) (Monitor Mode)
Temperature Predict Time	Oral: Approximately 4-6 seconds Adult Axillary (18 years and older): Approximately 12 - 15 seconds Pediatric Axillary (17 years and younger): Approximately 10 - 13 seconds Rectal: Approximately 10 - 13 seconds
Unit Dimensions	8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm)
Weight	12.6 oz (357 g)
Ambient Temperature Range	50.0 °F (10 °C) to 104 °F (40.0 °C)
Patient Temperature Range	80.0 °F (26.7 °C) to 110 °F (43.3 °C)
Power Source	(3) 1.5 Vdc AA batteries
Battery Life	Approximately 5,000 temperatures

PRODUCT AND ACCESSORY ORDERING INFORMATION

THERMOMETERS

01692-200: SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 4ft cord

01692-201: SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 4ft cord **01692-300**: SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord (must be locked to the wall)

01692-301: SureTemp Plus 692 Electronic Thermometer, One Per Room, Wall Mount, Rectal Probe and Rectal Probe Well with 9ft cord (must be locked to the wall)

01692-700: SureTemp Plus 692 Electronic Thermometer, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord, with Rolling Stand and Holder for Spare Probe & Well

ACCESSORIES

21326-0000: Wall Holder 21330-0000: Holder for Spare Probe and Well 06138-000: Calibration Key Assembly, 690/692 02892-000: Probe and Well Kit, 4ft, Rectal 02892-100: Probe and Well Kit, 9ft, Rectal 02893-000: Probe and Well Kit, 4ft, Oral 02893-100: Probe and Well Kit, 9ft, Oral 02891-0000: Probe Well, Blue 02891-1000: Probe Well, Red 08273-000: Thermometer Stand with Holder for Spare Probe and Well 21333-0000: Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well 98025-0000: Extended Warranty M692

REPLACEMENT PROBE COVERS

05031-101: Case of 1,000 Probe Covers

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

Baxter, Hillrom, SureTemp, and Welch Allyn are trademarks of Baxter International, Inc. or its