



AM-368 Super Duty Bariatric Treatment Table

- | | | | |
|--|---|---|--|
| <ul style="list-style-type: none"> • 3-Section Treatment Table • 36"W x 76"L | <ul style="list-style-type: none"> • Min. Height: 22" • Max Height: 38" • 1 3/4" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 36" with Contoured Face Opening • Center: 22" x 36" • Foot: 39" x 36" | <ul style="list-style-type: none"> • 800 lb. Capacity • 4 Total Lock Swivel Casters |
|--|---|---|--|



AM-F3 Treatment Table

- | | | | |
|--|---|---|---|
| <ul style="list-style-type: none"> • 3-Section Treatment Table • 25"W x 76"L | <ul style="list-style-type: none"> • Fixed Height: 30" • Shelf: 15.5" x 54" (upholstered) • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 25" with Contoured Face Opening • Center: 22" x 25" • Foot: 39" x 25" | <ul style="list-style-type: none"> • 400 lb. Capacity • Non-Skid Feet |
|--|---|---|---|