Avante Millennium 5

Hospital Stretcher

Available with easy-to-steer fifth wheel system

BENEFITS

- High-low adjustable mattress range is 13.5"
- Back-rest lifting range is 0 75° and controlled by oil hydraulics
- Knee-rest lifting range is 0 30° and controlled by adjustable crank
- Central braking system with steering function

- > Trendelenburg range is ±18°
- Two IV poles (one at the head and foot)
- 4" high density PU foam mattress covered with synthetic leather
- Collapsible side rails on both sides that can fold manually
- » 8" steering wheel central brake on four casters

FEATURES

- > 4" x 28.5" x 71.5" Pad
- > 8" Casters
- > 27" Wide Patient Surface
- Fold-Down Side Rails; Side Rails 14.5" x 39"
- > Stationary Foot Section
- > Four Brake/Steer Pedals
- > Dual Pedestal Hydraulics
- Dual Side-Mounted Foot Control Hydraulics
- Hands-Free Trendelenburg and Reverse Trendelenburg
- Steel bed frame with powder coating

- Mattress base made of ABS plastic board
- Bumpers made of durable ABS plastic and located at each corner
- Base cover made of ABS plastic with space for an oxygen bottle and emergency devices



SPECIFICATIONS



Weight Capacity: 550 lbs



Overall Length: 83"

Overall Width: 33"

Mattress Surface: 71.5" (L) x 28.5" (W)

Height Adjustment Range: 22.5" - 36.5" with 8" Casters



Trendelenburg/
Reverse Trendelenburg: ±18°

Backrest Adjustments: +90°



Casters:

8" with central locking one tracking 4 brake/steer pedals

I.V. Attachment Sockets: 2

The fifth wheel system saves power and allows for easy maneuverability

397-2019-10-29

