



Biological Monitoring

MAIL-IN AND IN-OFFICE SPORE TESTING SOLUTIONS

Helping dental offices cost effectively manage spore testing while maintaining CDC biological monitoring compliance.



BIOlogical 52 & 12 Mail-In Spore Test Service

AUTOMATIC TEST TRACKING WITH MINIMAL EFFORT

HealthFirst provides you with the most cost-effective way to perform your weekly autoclave spore testing, while providing 3rd party documentation of testing and test results. The included spore strips contain highly resistant dual-species bacterial spores providing the most reliable method of testing for an effective sterilization cycle.

Maximize Your Savings with Mail-In Solutions

Save money and stay compliant by using a mail-in solution for your weekly autoclave spore testing. Each mail-in test packet includes 1 test strip and 1 control strip, which satisfies CDC weekly biological monitoring recommendations. Test packets are compatible with steam, dry heat, chemical and ETO sterilizers.

OnTraq

OnTraq, HealthFirst's practice readiness application, reduces the risks, headaches and costs of biological monitoring regulatory compliance. It dynamically records spore test results for each sterilizer in your office. Test failures are immediately communicated via email and phone.

Rely on OnTraq's smart automation for dashboards, digital processing certificates, compliance reports and task management.



BIOlogical 52

52 Mail-In Test Envelopes
Item No. 1003600

BIOlogical 12

12 Mail-In Test Envelopes
Item No. 1003580



FREE SHIPPING! on select items
ordered on Auto-Replenishment*

HealthFirst.com | 800.331.1984

Confirm® 10 In-Office Spore Testing System

FINAL SPORE TEST RESULTS IN JUST 10 HOURS

Simply process the indicator in a normal sterilizer cycle. Activate the indicator, place it in the incubator, then see and record the results 10 hours later! If any bacterial spores survive the process, the indicator media will turn yellow, which indicates a failed test. Spore growth can be detected within 3-5 hours, giving you time to sterilize and re-test your equipment if issues arise.



Confirm 10 System Components:

Confirm 10

Steam Biological Indicators
25 Indicators / box
Item No. 1003610

Confirm

In-Office Monitoring Incubator
Incubator & manual
Item No. 1003590



FREE SHIPPING! on select items ordered on Auto-Replenishment*

Complete In-Office Testing System

Purchase a box of the Confirm 10 biological indicators (item 1003610) along with the re-usable incubator (item 1003590) for a complete in-office biological monitoring system for steam sterilizers. Use the OnTraq online application to log your test results for a complete reporting environment.

Providing you with infection control and dental waste management solutions



iBarrier

The Crosstex line of surface barrier film offers unique single-use disposable tablet computer sleeves.



Sharps/ Medical Waste Disposal

Disposal of sharps and medical waste utilizing cost-effective mailback solutions.



Pharmaceutical Waste

Disposal of unused and partially-used pharmaceutical waste, including dental carpules.

Your partner for practice readiness



HealthFirst.com | 800.331.1984