

Anesthesia and Respiratory Catalog | Airway Management

Gibeck® Trach-Vent®

The Trach-Vent HME and accessory products are designed specifically for tracheostomized or intubated, spontaneously breathing patients.

- Recommended for tidal volumes of 50 mL and above
- For use on both pediatric and adult patients
- Weighs only 4.0 g with 10 mL of dead space
- Trach-Vent holder conveniently attaches Trach-Vent to tracheostomy or endotracheal tube, allowing it to be removed and replaced easily during suctioning or change-out
- Oxy-Vent® snaps over Trach-Vent to allow delivery of supplemental oxygen



Trach-Vent® HME

41112 Item Code 41112

Individually packaged, 50 to a case.

Starter Kit

Item Code 41172

- 3 Trach-Vent HMEs
- 1 Oxy-Vent
- 1 Trach-Vent Holder

Individually packaged, 6 kits to a case.



41312

Trach-Vent®+ HME

Item Code 41312

- Compact, lightweight HME for spontaneously breathing patients
- Includes convenient, self-opening/closing suction port
- Integral oxygen port allows delivery of supplemental oxygen when needed
- Weighs only 6.5 g

Individually packaged, 50 to a case.

Trach-Vent® Heat & Moisture Exchangers

ITEM CODE	DESCRIPTION	TV RANGE (ML)	MOISTURE OUTPUT (MG H2O/L)	WEIGHT (G)	RESISTANCE (CM H2O)	DEAD SPACE (ML)	CASE QUANTITY
41112	TRACH-VENT	>50 mL	27, $Vt = 0.5 L$	4	0.25, 30 lpm	10	50
41312	TRACH-VENT+	>50 mL	27, $Vt = 0.5 L$	6.5	0.25, 30 lpm	10	50

© 2017 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions; contact customer service to confirm availability in your region.

Öãa[][¦óásæ) ÁÖða áæk\æ) ÁU|^@ÁK ÚVEKÖÖT SOTÁÚOESVQ R|ÁVæ) æ©ÁOGaæ) * ÁQQÓAÞ[ÁFJÁOEEÓ Rækæ ææÁEÁQÞÖUÞÒÙQOE Þ[ÁÜ^*ÁKASÖT ÒÞSÒÙÁÜQÁQOESŠÁGÆI €HFF€€Ì€