



REORDER NUMBERS:

10551 – WITHOUT HANDLE

10553 – WITH HANDLE

BARIATRIC FOOT STOOLS

500 LB. WEIGHT CAPACITY

The Dynarex Bari+Max Foot Stool is a high-quality, sturdy stool designed to safely provide additional height for larger individuals. Constructed with a heavy-duty steel frame and cross brace for added stability, this Foot Stool features reinforced rubber tips to help prevent slippage for optimal safety.

- Heavy-duty steel frame with cross brace provides optimal support
- Non-skid rubber platform for improved stability
- Tubular legs feature reinforced rubber tips to help prevent slippage
- Black powder coat finish for long-lasting protection
- 10553: Built-in handle provides additional stability and balance
- Limited Lifetime Warranty



10551



10553