



Polycarbonate Basket

Model 69576X1-P



Clear polycarbonate with rigid one-piece construction. It is available for all Novum bassinets.

Features:

- 3"x 5" card holder at head end
- Vent holes on bottom
- Dimensions: 26.5"L x 13"W 10"H

Bassinet Mattress

Model 7028X21



Novum mattresses conform to all federal regulations. The outside cover is stain&odor resistant, anti-bacterial & flame retardant. Foam material meets California #117 code for flame retardant material. All mattresses are latex-free. Additional sizes available.

Dimensions: 25.75"L x 12"W x 2"H