

SKIN FRICTION PROTECTORS







LiquiCell is a registered trademark and Shear Point is a trademark of Birchwood Laboratories LLC

Protects skin. Stays in Place.

LiquiCell[®] is a medical technology that helps prevent skin tissues from being stretched or torn. The ultrathin and nearly invisible liquid-filled membranes of LiquiCell[®] Shear Points, provide a synthetic, external bursa that allows the body to glide or float with an object both aiding in the reduction of shear and friction of the skin.

Unique Features:

- Self-adhesive Shear Points can be applied to the bridge of the nose, orbits of the eyes, ears, back of the heel or any other area at risk of skin breakdown due the presence of shear and friction
- Excellent for use with face masks, CPAP masks, Oxygen masks, tubing, glasses (eye/safety), podiatry (heel/toes) sports, braces, casts or any other article that contacts the skin
- Gentle skin-sensitive adhesive keeps LiquiCell Shear Points securely in place
- Transparency allows for easy monitoring of skin condition without removal
- Latex-free/Mercury-free

| - | SIZE | | QUANTITY | ITEM NO |
|---|--------|------|----------|---------|
| | Small | 56mm | 100 Ct | 94201 |
| | Medium | 70mm | 100 Ct | 94204 |



MADEIN

Birchwood Laboratories LLC 7900 Fuller Road, Eden Prairie, MN 55344-2195 800.328.6156 x7968 | Direct: 952.937.7968 medical@birchlabs.com | www.birchlabs.com