



## Birthing Mirror

*Model NBS-M1000N*



### Features:

- Heavy gauge steel for strength and durability
- Base with durable powder coat finish allowing for easy cleaning
- 5-leg, low-profile base is heavy weight and offers extra stability
- Extendable pole reaches a maximum of 70.5" in height
- Smooth rolling casters allow for accurate steering
- Easily adjustable mirror positioning

### Specifications:

- Mirror: 15" x 20.5"
- Height: 51" - 70.5"
- Base Diameter: 22"
- 3" Swivel Casters

### Base Color Options:

- Standard Black
- White
- Blue
- Red
- Yellow