



5-Year Parts, 1-Year Labor Warranty on cabinets  
 \*U.S./Canada only

**ENTHERMICS MEDICAL SYSTEMS**  
 An ISO 13485:2016 Certified Company  
 Tel: 262-251-8356 | 800-TO-B-WARM  
 generalinfo@enthermics.com  
 www.enthermics.com

EB1002 (02/2023)

### EC250

- Capacity: 4-6 blankets
- Interior size: 2.5 ft<sup>3</sup>
- Exterior size: 22" H; 18.5" W; 23.5" D
- Countertop model
- Stackable
- No fans or moving parts



### EC350

- Capacity: 10-12 blankets
- Interior size: 3.5 ft<sup>3</sup>
- Exterior size: 28" H; 18.5" W; 23.5" D
- Floor or countertop model
- Stackable
- No fans or moving parts



### EC400

- Capacity: 12-16 blankets
- Interior size: 4.0 ft<sup>3</sup>
- Exterior size: 22" H; 24" W; 25.75" D
- Stackable
- No fans or moving parts



### EC750

- Capacity: 20-25 blankets
- Interior size: 7.5 ft<sup>3</sup>
- Exterior size: 35.5" H; 24" W; 25.75" D
- Stackable
- No fans or moving parts



### EC1540

- Capacity: 40-50 blankets or 60 1-ltr irrigation bottles
- Interior size: 15.4 ft<sup>3</sup>, 7.7 ft<sup>3</sup> per chamber
- Exterior size: 72.75" H; 24" W; 32.5" D
- No fans or moving parts
- HCAI - OPM-0067-13



Shown with optional bumper

### EC2060

- Capacity: 60-70 blankets
- Interior size: 20.6 ft<sup>3</sup>
- Exterior size: 71" H; 32" W; 29.25" D
- HCAI - OPM-0067-13



## MULTI-PURPOSE WARMER

### EC550BL-18

- Total capacity:
  - TriMode chamber: 24 1-liter bottles, 20 1-liter bags or 14-18 blankets
- Exterior size: 25.12" H; 30.0" W; 21.12" D
- 18" chamber depth
- Also available as a built-in cabinet



Unit shown with fluid bottles, bags and blankets for illustration purposes only

### EC550BL-24

- Total capacity:
  - TriMode chamber: 30 1-liter bottles + 36 0.5-liter bottles, 48 1-liter bags or 16-20 blankets
- Exterior size: 25.12" H; 30.0" W; 28.5" D
- 24" chamber depth
- Also available as a built-in cabinet



Unit shown with optional basket - set up for fluids

## COMBINATION BLANKET & FLUID WARMERS (DUAL CAVITIES)

### EC1350BL

- Total capacity:
  - Blanket chamber: 24-30 blankets; 9.39 ft<sup>3</sup> interior size
  - TriMode chamber: 24 1-liter bottles, 20 1-liter bags or 9-12 blankets
- Exterior size: 75" H; 30.0" W; 21.25" D
- 18" chamber depth
- Also available as a built-in or pass-through cabinet
- OSHPD - OPM-0562-19



Shown as optional pass-through cabinet

### EC1850BL

- Total capacity:
  - Blanket chamber: 30-40 blankets; 13ft<sup>3</sup> interior size
  - TriMode chamber: 30 1-liter bottles, 24 1-liter bags or 10-14 blankets
- Exterior size: 75" H; 30" W; 28.5" D
- HCAI - OPM-0067-13
- \*height with casters



CONVENIENT CARE & COMFORTABLE PATIENTS

## STORAGE CABINET

### 1560-GS

- Interior size: 15.6 ft<sup>3</sup>
- Exterior size: 75" H; 30" W; 21.25" D
- Pairs with EC1350BL Combination Warmer



## BLANKET WARMERS

YOU CAN AFFORD THE BEST





ARE YOUR BLANKETS WARM ENOUGH?

# COMFORT SERIES

WARMSAFE ZONE HEAT – WARMS SAFELY, QUIETLY & EVENLY

## DC250

- Capacity: 4-6 blankets
- Interior size: 2.5 ft<sup>3</sup>
- Exterior size: 22" H; 18.5" W; 22" D
- Stackable
- No fans or moving parts



## DC350

- Capacity: 10-12 blankets
- Interior size: 3.5 ft<sup>3</sup>
- Exterior size: 28" H; 18.5" W; 22" D
- Stackable
- No fans or moving parts



## DC400

- Capacity: 12-16 blankets
- Interior size: 4.0 ft<sup>3</sup>
- Exterior size: 22" H; 24" W; 25" D
- Stackable
- No fans or moving parts



## DC750

- Capacity: 20-25 blankets
- Interior size: 7.5 ft<sup>3</sup>
- Exterior size: 35" H; 24" W; 24.3" D
- Stackable
- No fans or moving parts



## MIX & MATCH

Blankets, irrigation fluids and injection fluids need to be warmed to different temperatures. Match the right temperature to the right product in the right place at the right time. Only Entermics offers such a full line of stackable fluid and blanket warmers to create the warming configuration to fit your needs.



## STACKED BLANKET & FLUID WARMERS

### DC250L over DC350

- Capacity: 10-12 blankets and 12 1-liter bottles or 16 1-liter bags
- Exterior size: 50" H; 18.5" W; 22" D
- Additional combinations available



### EC400 over EC400L

- Capacity: 14-18 blankets and 28 1-liter bottles or 20 1-liter bags
- Exterior size: 44" H; 24" W; 24.3" D
- Additional combinations available



## MOBILE CARTS

### MS-400

- Mobile cart for small blanket and fluid warmers
- Fits the following models: DC250, DC250L, DC350, DC400, DC400L, EC250, EC250L, EC350, EC350L, ivNext®
- Exterior size: 36.5" H; 26" W; 26" D



Shown with DC250 Blanket Warmer

### MS-550

- Mobile cart for Multi-Purpose Warmers
- Fits the following model EC550BL-18" Depth, EC550BL-24" Depth
- Exterior size: 34" H; 32.5" W; 29.5" D



Shown with EC550BL 18" Depth Combination Warmer

## BLANKET WARMERS

WarmSafe™ zone heat – warms safely, quietly & evenly

Provide the warmest possible blankets and fluids with Comfort & Designer Series blanket warmers. These convenient, patient-friendly warmers are perfect for point-of-care placement.


All blanket warmers utilize WarmSafe zone heating to ensure every blanket is safely and evenly warmed. With no fans or moving parts on the tabletop blanket warmers, they are engineered for years of reliable operation.



DC250 with optional lock



EC750 with optional casters



Additional combinations available

\*Patent No: US 8,217,316; US 8,581,152