



www.boviemed.com

# **Suggested Procedural Uses for Bovie Cauteries**

### **High Temperature Fine Tip Cauteries**

	<b>AA01, H101</b> (pictured)	
Model		
Primary Use:	Pinpoint Haemostasis	Delicate, oculoplastic, dermal, and minor subcutaneous procedures.
By Specialty:	Cardiology	Carotid Endarterectomy, Valve Replacement, Femoral Popliteal Artery Bypass, Patent Ductus Arterlosus, Saphenous Vein Ligation, Arteriovenous Shunt, Fistula Formation, Pacemaker Implantation.
	Neurosurgery	Craniotomy, Tumor Removal, Acoustic Neuroma, Peripheral Nerve Repair.
	OB/GYN	Abdominal Hysterectomy, Mastectomy, Vulvectomy, Perineal Lacerations, Vaginoplasty, Salpingectomy, Epislotomy Repair, Cesarean Section, Repair of Third Degree, Circumcision.
	Ophthalmology	Strabismus Surgery, Tumor and Cyst Removal, Enucleation, Keratectomy, Blepharoplasty, Dacryocystorhinostomy.
	Orthopaedics	Total Hip, Total Knee, Arthroplasty, Arthrotomy, Open Reduction of Fractures, Amputation of Extremities, Shoulder Fusion, Spinal Fusion, Bone Graft Procedures.
	Otolaryngology	Thyroidectomy, Tumor Removal, Sinus Surgery, Tracheostomy, Bepharoplasty, Carbuncle Surgery, Parotid Gland Excision, Radical Neck Dissection, Rhinoplasty, Nasal Septal Surgery.
	Plastic Surgery	Blepharoplasty, Mammoplasty, Rhinoplasty, Septonasal Reconstruction, Submucous Resection of Septum, Cleft Palate, Rhytidectomy, Abdominal Lipectomy, Tendon Repair, Scar Revision, Lacerations, Hair Transplants, Burn Procedures
	Proctology	Hemorrhoidectomy





www.boviemed.com

## **Suggested Procedural Uses for Bovie Cauteries**

#### **Low Temperature Fine Tips**

Model	AA00, H100 (pictured)	
Primary Use:	Pinpoint Haemostasis	Delicate, oculoplastic, dermal, and minor subcutaneous procedures.
By Specialty:	Neurosurgery	Craniotomy, Peripheral Nerve Repair, Surgery on Aneurisms
	Ophthalmology	Cataract Extraction, Corneal Transplant, Repair of Burns and Lacerations, Strabismus Surgery, Tumor and Cyst Removal, Enucleation, Retinal Surgery, Keratectomy, Dacryocystorhinostomy
	Plastic Surgery	Tendon Repair

#### **High Temperature Loop Tips**

Model	AA03, H103 (pictured)	
Primary Use:	Haemostasis	Hemostasis of small bleeding areas, incision and excision
By Specialty:	Neurosurgery	Tumor Removal
	OB/GYN	Incision and Drainage of Abcess, Removal of Vulvular Condyloma, Repair of Third Degree.
	Ophthalmology	Strabismus Surgery, Enucleation, Keratectomy.
	Orthopaedics	Total Knee, Arthroplasty, Arthrotomy, Open Reduction of Fractures, Amputation of Extremities, Bonlonectomy.
	Otolaryngology	Thyroidectomy, Tumor Removal, Carbuncle Surgery, Parotid Gland Excision, Division of Scalenus Anticus Muscle.
	Plastic Surgery	Cleft Palate, Scar Revision, Lacerations, Hair Transplants, Burn Procedures.
	Proctology	Surgery on Rectal Condyloma.
	Urology	Proud Flesh and Veneral Wart Removal



# **Suggested Procedural Uses for Bovie Cauteries**

### **High Temperature 2" Extended Loop Tips**

	<b>AA05, H105</b> (pictor	ured)
Model		4)
Primary Use:	Control of small bleeding areas in recessed or limited access places, skin lesion removal.	
	Cardiology	Valve Replacement, Ventricular Septal Defect (VSD), Aortic Aneurism, Tetralogy of Fallot, Femoral Popliteal Artery Bypass, Patent Ductus Arterlosus, Saphenous Vein Ligation, Arteriovenous Shunt, Fistula Formation.
	Neurosurgery	Acoustic Neuroma, Laminectomy.
	OB/GYN	Cervical Polyp, Trachelorrhaphy, Abdominal Hysterectomy, Mastectomy, Tubal Ligation, Anterior and Posterior Repair of Vaginal Walls, Incisions and Drainage of Abcess, Removal of Vulvular Condyloma, Vaginoplasty, Salpingectomy, Episiotomy Repair, Cesarean Section, Repair of Third Degree, Laceration, Repair of Cystocele, Repair of Rectocele.
	Ophthalmology	Blepharoplasty, Dacryocystorhinostomy
By Specialty:	Orthopaedics	Total Hip, Total Knee, Arthroplasty, Arthrotomy, Open Reduction of Fractures, Amputation of Extremities, Shoulder Fusion, Spinal Fusion, Bone Graft Procedures
	Otolaryngology	Thyroidectomy, Nasal Polyps, Tumor Removal, Sinus Surgery, Control of Eplstaxis, Tracheostomy, Bepharoplasty, Parotid Gland Excision, Radical Neck Dissection, Rhinoplasty, Nasal Septal Surgery, Foreign Body Removal.
	Plastic Surgery	Blepharoplasty, Mammoplasty, Rhinoplasty, Septonasal Reconstruction, Submucous Resection of Septum, Cleft Palate, Rhytidectomy, Abdominal Lipectomy, Nasal Polyps, Burn Procedures.
	Proctology	Colon Resection, Rectum Resection, Rectal Condyloma
	Urology	Hydrolectomy



# **Suggested Procedural Uses for Bovie Cauteries**


## **High Temperature 5" Extended Loop Tips**

del		
mary Use:	Control of small bleeding areas, incision and excision.  5" tip for places requiring extended reach.	
	Cardiology	Valve Replacement, Ventricular Septal Defect (VSD), Aortic Aneurism, Tetralogy of Fallot, Femoral Popliteal Artery Bypass, Patent Ductus Arterlosus, Saphenous Vein Ligation, Arteriovenous Shunt, Fistula Formation.
	Neurosurgery	Acoustic Neuroma, Laminectomy.
	OB/GYN	Cervical Polyp, Trachelorrhaphy, Abdominal Hysterectomy, Mastectomy, Tubal Ligation, Anterior and Posterior Repair of Vaginal Walls, Incisions and Drainage of Abcess, Removal of Vulvular Condyloma, Vaginoplasty, Salpingectomy, Episiotomy Repair, Cesarean Section, Repair of Third Degree, Laceration, Repair of Cystocele, Repair of Rectocele
	Ophthalmology	Blepharoplasty, Dacryocystorhinostomy
Specialty:	Orthopaedics	Total Hip, Total Knee, Arthroplasty, Arthrotomy, Open Reduction of Fractures, Amputation of Extremities, Shoulder Fusion, Spinal Fusion, Bone Graft Procedures.
	Otolaryngology	Thyroidectomy, Nasal Polyps, Tumor Removal, Sinus Surgery, Control of Eplstaxis, Tracheostomy, Blepharoplasty, Parotid Gland Excision, Radical Neck Dissection, Rhinoplasty, Nasal Septal Surgery, Foreign Body Removal.
	Plastic Surgery	Blepharoplasty, Mammoplasty, Rhinoplasty, Septonasal Reconstruction, Submucous Resection of Septum, Cleft Palate, Rhytidectomy, Abdominal Lipectomy, Nasal Polyps, Burn Procedures.
	Proctology	Colon Resection, Rectum Resection, Rectal Condyloma.
	Urology	Hydrolectomy



## **Suggested Procedural Uses for Bovie Cauteries**

#### **High Temperature 5" Extended Fine Tips**



#### **High Temperature Vasectomy & Micro Vasectomy Tips**

Model	<b>AA21, AA29, H121</b> (pictured)	
Primary Use:	Destruction of the mucosa inside the vas deferens, sculpting of woven skin grafts.	
By Specialty:	Urology	Intravessicle Removal of Stones, Hydroelectomy, Vasectomy.