accuc@ld

ARG49ML



83.75" x 55.25" x 31" (H x W x D)

Performance series Pharma-Lab 49 cu.ft. all-refrigerator in stainless steel with glass doors

Highlights:

Intelligent microprocessor digital temperature controller

Improved temperature accuracy with a separate temperature probe for continuous product readout

Additional sensor reads interior air temperature for superior temperature control

Product Features:

Purpose-built design	Designed and engineered to store vaccines and other sensitive materials under stable temperature conditions	
Glass doors	Full display of the interior, with an anti-condensate door frame heater to ensure a clear view	
Adjustable temperature	Interior is pre-set to 5°C, with an adjustable control range from 2 to 10°C	
Temperature alarm	Visual and audible alarms alert if the interior temperature rises or falls out of set range	
Complete stainless steel construction	Interior and exterior are entirely constructed from stainless steel for lasting durability and easy cleanup	
Forced air cooling	Top-mounted commercial fan ensures rapid temperature recovery	
Automatic defrost with adjustable cycle defrost	Adjustable cycle defrost operation ensures the unit avoid temperature spikes and ice buildup	
Self-closing doors	Two doors close to save energy but will remain open when the door is opened at a 90° angle for easier loading	
Open door alarm	Alarm will beep 3 times if the door is left open for 30 seconds, 5 times at 60 seconds, and will sound continuously if the door is open for more than 5 minutes	



accuc@ld

ARG49ML Specifications:

Overview		
Height of Cabinet	83.75" (213 cm)	
Width	55.25" (140 cm)	
Depth	31.0" (79 cm)	
Depth with Handle	32.75" (83 cm)	
Depth with door at 90°	56.5" (144 cm)	
Capacity	49.0 cu.ft. (1388 L)	
Defrost Type	Automatic	
Door	Glass	
Cabinet	Stainless Steel	
US Electrical Safety	UL	
Canadian Electrical Safety	ULC	
Sanitation	ETL-S	
Amps	7.5	
Voltage/Frequency	115 V AC/60 Hz	
Weight	410.0 lbs. (186 kg)	
Shipping Weight	516.0 lbs. (234 kg)	
Parts & Labor Warranty	2 Years	
Compressor Warranty	5 Years	
UPC	761101062341	
Refrigerator Features		
Door Swing	French	
Shelf Type	Wire	
Shelf Qty	8	
Adjustable Shelves	Yes	
Thermostat Type	Digital	
Fan Type	Interior	
Refrigerant Type	R290a	
Refrigerant Amount	3.2oz.	
High Side PSI	290.0	
Low Side PSI	155.0	
Wheel Qty	4	
Interior Light	Yes	
Temperature Range	2 to 10°C	
Dimensions		
Interior Height	51.5" (131 cm)	
Interior Width	50.5" (128 cm)	

