
# PRODUCT DATA SHEET




# ResQPod

THE RESQPOD ENABLES CARERS TO RAPIDLY AND SAFELY EVACUATE PEOPLE WITH MOBILITY PROBLEMS IN AN EMERGENCY.

The ResQpod is an evacuation board designed to help evacuate non-ambulant people or people who may not be able to walk down stairs in an emergency. It can be used for both children and adults.

The ResQpod can be used in office buildings, theatres, libraries, hospitals, schools and many other facilities.



#### **Features**

#### 1 Long pulling straps

Located at both head and feet to allow for multiple users

#### 2 Secure body protection

Completely encloses the upper body (with the option of having arms out)

#### 3 Heavy duty hook and loop

For increased pssenger safety

### 4 Impact absorbing pillow

For added head protection when transporting down stairs or across rough terrain

#### 5 Glideable plastic board

Corrugated plastic board allows for excellent heat insulation preventing injury, plus removing the necessity to unroll every few months to prevent the plastic keeping its shape

# 6 Adjustable Straps

Curls the ResQMat around the feet, allowing it to fit all sizes

#### 7 Padded mini mattress

Situated under the body for added comfort during transportation

#### **Includes**

- 440lbs Weight capacity
- 2 person operation
- Can be moved easily over virtually any surface
- Narrow width and flexibility
- Wipe clean storage bag
- Can be used down stairs

#### **Benefits**

- No installation necessary, ready to use instantly
- No belay system required
- No dangerous ropes to negotiate during evacuation

Dimensions	pac	ked

Height	20"
Width	24"
Depth	5"
Weight	4.8lbs

## Dimensions unpacked

Height	89"
Width	51"
Depth	3"
Weight Capacity	440lbs

# VIDEO LINK: evac-chair.com/evacuation-sleds/resqpod/

For over 35 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®.

Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11, 2 YEAR MANUFACTURERS WARRANTY.

For further information please contact:

